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The device  
harmonizes  
and corrects the natural  
electromagnetic emissions  
of natural bio systems

DeVita Ritm  
Model Base

## Introduction

For the first time the thought of energy presence in a human body and of health dependence on its circulation was expressed by ancient Esculapians of China, India, and Egypt. Every human organ emits certain resonance frequencies. Under the influence of unfavorable external factors, modern life full of stresses, bad environment, its frequencies deviate from their natural parameters, and a risk of disease can occur.

The «DeVita Ritm Model Base» electromagnetic device is a unique development in the field of supporting good health without medications. “DeVita Ritm Model Base” softly and safely influences the human body without disturbing its energetic balance either in the whole or locally.

The most important advantage of “DeVita Ritm Model Base” is the availability of automatic programs for health support by sessions. In addition, the device has a feature of programming by any of the available prophylactic programs with «DeInfo». This allows you to choose an individual prophylactic course for everyone.

## List of programs

- |  |                                 |
|--|---------------------------------|
| 1. 7 Chakras                             | 18. Intestine without problems  |
| 2. Active protection                     | 19. Joints without inflammation |
| 3. Anti-Cold                             | 20. Joints without pain         |
| 4. Anti-Stress                           | 21. Perfect eyesight            |
| 5. Blood circulation (Regulation)        | 22. Perfect hair & nails        |
| 6. Cardiac rhythm normalization          | 23. Perfect mood                |
| 7. Clean kidneys                         | 24. Perfect skin                |
| 8. Clean lymph                           | 25. Recovery from injuries      |
| 9. Detoxification                        | 26. Sexuality, Female           |
| 10. Digestion                            | 27. Sexuality, Male             |
| 11. Digestion (Metabolism) normalization | 28. Sports activities. Workout. |
| 12. Free breath                          | 29. Warm-up                     |
| 13. Health, Female                       | 30. Without allergens           |
| 14. Health, Male                         |                                 |
| 15. Healthy heart                        | <b>Automatic Programs:</b>      |
| 16. Healthy sleep                        | 1. Healthy day                  |
| 17. Hormones normalization               | 2. Complete cleansing           |

**Caution!** The device has contraindications for use by women of the first trimester of pregnancy and by people who underwent organ transplantation. It is also recommended not to use the device during the first two months after myocardial infarction. In case of severe heart diseases it is necessary not to place the device closer than 0.5 m from the body.

## 1. 7 Chakras

The human being is a complex energy structure which is based on seven Chakras located along the spine. Each chakra is associated with a certain organ, thus, when chakras are damaged, the function of their respective organs is violated. The chakras accumulate energy. Negative energy clogs up the chakras and has a detrimental effect on the whole body.

The program recovers lost energy, correctly allocates it and thereby restores the organs and the systems.

**Use:** The program is used during the daytime.

## 2. Active protection

The human body exists in a constantly changing environment that it has to adapt to. It is negatively affected by changes in temperature and atmospheric pressure, electromagnetic radiation of household appliances and office equipment, Along with frequent stress. The body is also affected by physical and mental load and ever-changing work pattern.

The program normalizes the nervous and endocrine systems, improves attention and capacity for work, mitigates fatigue and irritability and keeps you on your toes. It protects your organism against pathogenic radiation.

**Use:** Daily, in the morning.

## 3. Anti-Cold

The first signs of a cold are a stuffy nose and sore throat. We often pay no attention to these symptoms and carry on whilst ill. It is fraught with complications, the most common of which are upper respiratory tract infections such as: antritis, sinusitis, bronchitis and otitis.

The program activates the immune system, increases resistance, reduces the likelihood of complications and swelling, ensures better discharge of mucus from the sinuses and helps recovery in a short time.

**Use:** At any time of the day; during the off-season or cold epidemics use daily. It relieves allergies.

## 4. Anti-Stress

Almost every day we experience stress: at work, on the road, at home. Psyche usually neutralizes the effects of stress, but with constant negative stratification or abrupt change in the rhythm of life, the body needs aid.

The program will help in a difficult situation, when psychological or physical stress occurs. It normalizes the nervous system and ensures rapid falling asleep and sound sleep. In case of strong nervous excitement, use it in conjunction with the program "Healthy sleep".

**Use:** Before bedtime. The device is placed under the pillow. If sleep does not come at the end of the program, it is necessary to repeat the session.

**Caution!** The mode is not recommended for use while driving.

## 5. Blood circulation (Regulation)

Sedentary lifestyle, bad habits, poor diet and stress have a negative impact on the state of the circulatory system. The program recovers local and general blood circulation and regulation of oxygen consumption, which improves blood flow to all the organs.

**Use:** At any time, whenever necessary.

**Caution!** It is not recommended to use the program within 2 months after myocardial infarction.

## 6. Cardiac rhythm normalization

When the heart beats as it should, we do not notice its work. But if it starts to stray from its rhythm it brings huge discomfort. Arrhythmia can occur in healthy people for various reasons, such as a reaction to tight clothing, insect bites, certain medications, stress, abundant food, problems with stool, etc.

The program restores heart rhythm, improves blood flow to the heart and normalizes the nervous system.

**Use:** The program is used once a day, when the violation of the heart rate takes place.

## 7. Clean kidneys

The functions performed by kidneys are diverse. Their main function is the regulation of water and electrolyte metabolism in the body to maintain a constant volume and pressure of blood and other fluids. Kidneys also regulate the acid-base balance in the body. They withdraw the products of nitrogen metabolism: urea, creatinine, uric acid, and others.

Eating and metabolic disorders, the excessive consumption of meat, lack of vitamin A, hypothermia and infectious contamination – cause pollution of the kidneys and the formation of sand and stones.

The program ensures active removal of toxins from the body, regulates the work of the kidneys, and prevents the formation of kidney stones.

**Use:** The program can be used at any time of the day, optimally, from 5 p.m. to 7 p.m. In case of infectious or toxic contamination of the body, repeat the program in

30 minutes. For the normal function of the kidneys, it is necessary to drink clean water – 1-1.5 liters per day.

**Caution!** During this program there may be an increased urine output.

## 8. Clean lymph

The lymphatic system is the guardian of the purity of the body around-the-clock. It plays an important role in the removal of metabolic products, as well as in the maintenance of immunity. The lymph is cleaned naturally when it passes through lymph nodes, but sometimes this is not enough. In cases where dysfunction of the lymphatic system occurs, up to 83% of harmful substances may accumulate in the intercellular space. This increases the stress to all the organs of elimination: liver, intestine and kidneys.

The program ensures the liquefaction of lymph and the removal of toxins; it removes excess water from the body, improves cell metabolism, normalizes blood pressure and reinforces the immune system.

**Use:** The program can be used at any time of the day. During the program it is recommended to massage the area with gentle strokes.

**Caution!** For the most efficient detoxification, drink clean water – 1-1.5 liters per day.

## 9. Detoxification

Throughout our lives toxins accumulate in our bodies. They are accumulated in both intercellular substance and the cells of organs. Toxins disrupt metabolic processes and create a favourable environment for the breeding of pathogens.

The program activates the work of the excretory systems, helps to cope with the withdrawal of toxins from the body and regulates the liver and gall bladder functions.

**Use:** At any time of the day. Be sure to drink clean water – 1-1.5 litres per day.

## 10. Digestion

Factors, such as dietary violation and shift to unfamiliar food products, have an increased load on the gastrointestinal tract. This causes indigestion, flatulence, discomfort and pain.

The program affects favourably on the gastrointestinal tract. It helps the body assimilate food more easily, decreasing the risk of eating disorders and reducing the feeling of heaviness and discomfort. It relieves heartburn and bowel spasms.

**Use:** The program can be used at any time of the day, optimally, from 7 p.m. to 9 p.m. If necessary, re-use in 30 minutes.

## 11. Digestion (Metabolism) normalization

Metabolism is a set of chemical reactions that ensure the normal functioning of the organism. The main factors of metabolic disorders are: poor diet, use of certain

medications, stress and sedentary lifestyle. Metabolic disorders can cause such diseases as diabetes, gout, obesity, allergies, disorders of the cardiovascular system, growth and developmental delay in children.

The program is aimed at a balanced function of thyroid and pancreas, proper assimilation of potassium, calcium, magnesium, sodium and alignment of acid-base balance.

**Use:** 1-2 times a day, if necessary.

## 12. Free breath

Bronchospasm (cough) is a protective reaction of the organism to a stimulus and respiratory tract infections; it may also be a consequence of the effects of certain medicines on the body. Bronchospasm is common among people with a tendency to allergies and those suffering from bronchial asthma.

The program restores the function of the bronchopulmonary system, normalizes lung function, regulates blood circulation and oxygen uptake and relieves spasms.

**Use:** Use the program daily, if necessary, at any time of the day. If coughs and spasms are strong – use it up to 2-3 times in a row.

## 13. Health, Female

Bad environment, bad habits and sunlight adversely affect the female body, accelerate the ageing process and cause health problems.

The program covers a wide range of female health' problems and is useful at any age. Decreases premenstrual tension and pain, and adjusts menopausal disorders.

**Use:** For the purpose of prevention, the program is used once a week, at any time. During puberty, girls are recommended to use it 1-2 times a week. For women, experiencing pain and discomfort during menstruation, it is recommended to use the program every day during this period. In the correction of problems in the sexual sphere and/or in violation of the organs of the female reproductive system, for an effective impact use it on a daily basis in conjunction with the program "Sexuality, Female".

**Caution!** It is not recommended to use the program during pregnancy or if you are trying to conceive.

## 14. Health, Male

Male health is not only about healthy potency, but also about the correct functioning of other internal organs. For men, it is necessary to take care of their health from a young age.

The program is intended to harmonize the male urogenital system. It regulates urination process and the function of male sexual glands, including the prostate. The older a man, the more useful the program.

**Use:** Once a week, optimally, from 9 p.m. to 11 p.m. If you have problems with potency, and/or prostate, use it on a daily basis in conjunction with the program "Sexuality, Male".

## 15. Healthy heart

The supply of organs with oxygen and nutrients depends on the heart function working effectively, to support the whole organism. Failures in the activity of the heart are caused by unhealthy lifestyle, stress and physical overload.

The program adjusts the integrated circulation and blood pressure, provides normal blood supply to the heart muscle, and ensures self-regulation of the heart and its stable function.

**Use:** At any time of the day, once or twice, if necessary (optimally, from 11 a.m. to 1 p.m.).

**Caution!** It is not recommended to use the program within 2 months after myocardial infarction.

## 16. Healthy sleep

Healthy deep sleep is an important part of everyone's life. During sleep, the body restores energy spent during the day. However, any psycho-emotional upheaval may serve as an impetus for the development of sleep disorders. This in turn provokes the development of both physiological and mental disorders. Irritability, fatigue, inability to concentrate, are the very first signs of a sleep disorder. The program helps to reduce excitability of the nervous system, ensures the rapid falling asleep and normalizes the sleep phase.

**Use:** Daily at night, re-use of the program up to 2-3 times is possible.

## 17. Hormones normalization

The endocrine system regulates the function of organs with the help of hormones and provides a constant internal environment (homoeostasis). The endocrine system is involved in the body's development, its growth and reproductive function. It also regulates the mental state and emotional reactions.

The program normalizes the function of the system "hypothalamus – hypophysis – adrenal glands – reproductive glands" and regulates excretory processes in the body.

**Use:** The program can be used at any time of the day.

## 18. Intestine without problems

Who among us has never experienced such symptoms as distension, nausea, intestinal problems, loss of appetite, abdominal pain? All of them are signs of disorders in the intestine function. These symptoms are usually not revealed in early stages, but when the disease has become chronic. The program reduces pain,

eliminates nausea and distension, normalizes stool, as well as fights with a common modern problem that is dysbiosis.

**Use:** Daily, if necessary; depending on the state – 1-2 times per day.

## 19. Joints without inflammation

Every third inhabitant of our planet experiences pain or crunching in the joints and spine. The reasons for this problem could be infections, traumas, allergies, overweight, professional sports, poor diet, etc. The program improves the blood supply to the injured area and reduces inflammation, swelling and pain. It has revealed itself well as a supporting and restoring means during the “active garden work”, as well as in case of sprains and strong muscle fatigue. It has a favorable effect on the nerve endings of the joints and is used in case of pinching and inflammation of nervous tissue.

**Use:** Daily, 2 weeks prior to the active work in the garden and sports activities during a month, then – 1-2 times a week till the end of the season. For people with a periodic worsening of back and joints, in order to maintain the health – use it 2 times a year, daily, during a month. In case of exacerbation, combine with the program “Joints without pain”.

## 20. Joints without pain

Pain in the joints or spine can have several causes. Most often the result of pinched or inflamed nerves that occur due to a sedentary lifestyle, bruises, falls, overweight, herniated discs of the spine, hypothermia, active physical exercises, destruction of cartilaginous tissue and inflammation of the joint capsule. The program helps to reduce pain, improve blood microcirculation and normalize the muscle tone, ligaments and cartilage, and activates metabolic processes in damaged tissues.

**Use:** During exacerbations use the mode in addition to primary therapy, or in combination with the program “Joints without inflammation”. Every day, at any time of the day for 2 weeks, take a one day break, then repeat the course until you achieve the desired results.

## 21. Perfect eyesight

The constant strain on the eyes from reading or computer work, as well as the change of the conditions and the rhythm of life cause fatigue, eye irritation and a decrease of vision.

The program reduces the strain on the eyes, relieves fatigue and pain and ensures the normalization of vision.

**Use:** Before bedtime each day for 7 days. Then, one day of rest. Repeat the sessions until the desired result. Do not strain the eyes after the end of the program. For the best result, combine with exercises for the ciliary muscle.

## 22. Perfect hair & nails

Hair loss, dandruff, brittle nails are indicators of body function disorders, or adverse environmental effects. The program ensures the improvement of blood microcirculation, activates cell metabolism and cell regeneration. This leads to the strengthening of the hair follicles and nails, stimulating hair growth and an increase of protective functions of the skin.

**Use:** Daily, until you get the desired result.

## 23. Perfect mood

Depression, anxiety, fear... These are interrelated concepts. Anxiety and fear are signs of depression. This condition very often leads to psychosomatic illnesses, emotional exhaustion and dissatisfaction. Fears are strong factors that prevent a person from developing his/her capacities and achieving success. "Perfect mood" program stops depression, reduces anxiety and normalizes sleep, it improves mood and mental alertness, as well.

**Use:** Daily, once during the day.

## 24. Perfect skin

Skin problems occur due to the dysfunction of the internal organs. In this case, the skin takes on a cleaning function. Acne, vascular defects, pigmentation and other similar conditions appear on the face. These are faced by teenagers experiencing hormonal changes and persons suffering from disorders of the endocrine and digestive systems.

The program normalizes the skin glands, strengthens the skin barrier function and improves metabolic processes in cells.

**Use:** 1-2 times a day until you achieve the desired result.

## 25. Recovery from injuries

Bruises, sprains and fractures can happen to both sportsmen/women and ordinary people. In conjunction with medical aid, the program speeds up the healing process, reduces pain and swelling, eliminates bleeding and restores elasticity of muscles.

**Use:** Daily, 2-3 times a day until you achieve the desired result. If necessary, combine with the programs "Joints without inflammation" and/or "Joints without pain".

## 26. Sexuality, Female

Sexual desire is an indicator of various female states. If a woman is no longer interested in sexual activity, a disease can be assumed to have arisen. In most cases, a decreased sex drive develops due to a combination of several factors: physical, mental and microsocial ones.

The program will help in the absence or diminished sex drive in women; it will lessen the sense of emotional stress and mental discomfort, and enhance sexual desire.

**Use:** Daily, once a day.

## 27. Sexuality, Male

Manifestation of male drive depends on age, physical activity, nervous and mental state, external stimuli, etc. Stresses, problems at work and various complexes lead to psychological impotence in men. As a result, they become wary to the emergence and "quality" of erection, they compare and analyze. And the more a man prepares for this "event", the less his chances of success are.

The program is aimed at psychological adjustment and regulation of work of the 6th Chakra, it reduces anxiety and strengthens confidence in the favorable completion of intimacy.

**Use:** Daily, once a day.

## 28. Sports activities. Workout.

Successful and active workouts often lead to fatigue and pain in the muscles, which make it impossible to achieve greater results.

The program increases muscle strength, relieves muscle fatigue, strengthens the muscles, improves the supply of muscles with blood and oxygen, increases endurance, allows achieving maximum results with minimum effort of the body.

**Use:** During the workout.

**Caution!** It is not doping.

## 29. Warm-up

Gym-goers do not often have enough energy and motivation, especially after the working day. "Warm-up" program energizes, increases the intensity of training, it helps to overcome fatigue, improves metabolic processes in the muscles and makes the workout more effective.

**Use:** Before every gym session.

## 30. Without allergens

Allergens are the substances to which the body has an increased sensitivity due to the impaired immune response mechanism. They are all around us: outdoors, at home, at work. Common allergens are pollen, dust, animal dander, medications, household chemicals and others. The body can react to these substances with such symptoms as watery eyes, swelling of the mucous membranes, skin rash, etc. The program facilitates the allergy symptoms and normalizes the function of the immune system.

**Use:** On a daily basis, if necessary. For those who have allergies to flowering plants, it is recommended to start applying the program one month before the flowering.

## Automatic Programs

### 1. Healthy day

The complex "Healthy Day" improves the blood supply to organs and tissues, activates the excretory system, which, in turn, improves metabolism and correct distribution of energy balance, relieves nervous tension and strengthens immunity.

Includes the following programs: Active protection, Blood circulation (Regulation), Anti-Stress, Detoxification, Digestion, Hormones normalization, 7 Chakras.

**Use:** Daily in the morning, during a month. The use of the complex is most effective during colds, high physical, mental and psychological stress, the period of rehabilitation after illness or surgery.

### 2. Complete cleansing

We are constantly exposed to toxins. Some of them are produced by our body, others are received with food, water and air, and some of them are produced by pathogens that enter our body. The human body is endowed with a self-cleaning system by nature, but with modern ecology it is often unable to cope with its task. Slagging does not only worsen the general condition, but also paralyses the work of the organs and systems.

The complex cleans the system and organs from waste and toxins, regulates the function of the excretory organs, and improves the immune system.

Includes the following programs: Active protection, Blood circulation (Regulation), Intestine without problems, Clean kidneys, Clean lymph, Hormones normalization.

**Use:** Daily (optimally, from 8 a.m. to 11 p.m.) for 3-4 weeks, a two-week break, repeat the cycle again twice. For the purpose of prevention, use every six months, every day during a month.

**Caution!** While using the complex programs, sweating, odour, discolouration of urine and other manifestations of toxic withdrawal from the body are possible. Be sure to drink clean slightly alkaline carbonated water, at least 1.5 litres/day.

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DETA-ELIS-EUROPA GmbH